

WONDER BULK FLOW

THE SELF ADAPTATIVE FLOWABLE BULK FILL COMPOSITE

Several clinical situations require dentists to work swiftly against time, for instance;

Class II restorations with **deep gingival seat**

Class V cavities **close to gingival sulcus**

Patients with **excess crevicular fluid or saliva**

Such clinical scenarios do not usually give adequate time for proper instrumentation and manipulation.

Wonder Bulk Flow is an advanced self-levelling flowable composite designed to deliver a direct and quick application. With a 4mm depth of cure, it helps you complete complex restorations with ease and speed.



Packing Best in Class filler load (73%w/w) **Wonder Bulk Flow** offers amazing strength and does not require a capping layer on top. Just place the self-levelling resin directly from the syringe and watch it self-adapt to the cavity walls ensuring void-free restorations.

- Unique **Self-levelling** feature ensures close adaptation without much instrumentation and manipulation
- **Quick, easy, stress-free bulk placement** – helps minimize voids delivering more durable monobloc restorations



Best in class filler load of 73% ensures high compressive strength and abrasion resistance

No capping layer required

- **One-time cure of up-to 4mm** saves time
- Close adaptation to the cavity walls ensures **better marginal integrity** of the restoration.
- Easy to place in **Geriatric and Paediatric patients**

Equipped with special **symmetric refractive fillers**, no need to compromise on shade matching while using a bulk composite. Moreover, these fillers intelligently adapt to the optical density of surrounding tooth structure, hence achieving an appropriate shade match with a single universal syringe.

SPECIALLY INDICATED FOR



UP TO 4 mm
SINGLE STEP CURE

Bulk Posterior Restorations






Tunnel Restoration



Class V Abfractions

ORDERING INFORMATION

-  2g Syringe x 1
-  Applicator Tips x 5
-  1 Shade – Universal

[Click here](#) to find more on what other dentists are talking about Wizdent Products

 fb.com/fineartofdentistry

 www.wizdent.com