



# BONDING AGENT

# **SELECT BOND TE**

## THE INTUITIVE 5<sup>TH</sup> GEN DENTAL ADHESIVE

Post-Operative Hypersensitivity in composite bonding has challenged clinicians for many years. Among other factors impacting bond performance, passive evaporation of the solvent (ethanol / acetone) from the bottle has shown to profoundly impact bond viscosity & bond penetration. This results in a gap between the etched tubule and the resin tag; and often causes post-operative hypersensitivity.

**Select Bond TE's** all adhesive formulation prevents solvent evaporation and virtually eliminates post-operative hypersensitivity.





Select Bond TE's unique dual primer formulation enhances its wettability and helps the bond penetrate deeper into the dentinal tubules. This results in better resin tag formation and ensures stronger bonds to both dentin & enamel.

Unique **All adhesive No Solvent** formulation

Virtually **eliminates post-operative hypersensitivity**

Doesn't require air thinning -  
**Less technique sensitive**

**Saves time** - Just paint,  
wait and cure

Improves bond integrity by  
**minimizing surface contamination**

100% adhesive formulation guarantees  
**more restorations per bottle**

**Lasting bonds** – the trusted legacy of PIDILITE

Its unique solvent-free formulation eliminates the need for air drying/thinning making it fast and easy to use – just paint, wait and cure.



## SPECIALLY INDICATED FOR

- Adhesive for direct composite restorations
- Adhesive for bonding to composite and amalgams
- Adhesive for direct and indirect inlays, onlays post and veneers

## APPLICATION

Dispense directly on applicator tip

No air thinning required. wait 20 sec



Apply single coat



Light cure for 20 sec

## ORDERING INFORMATION

 5 ml Bottle

[Click here](#) to find more on what other dentists are talking about Wizdent Products

 [fb.com/fineartofdentistry](https://www.facebook.com/fineartofdentistry)

 [www.wizdent.com](http://www.wizdent.com)